

MapQ Tutorial

Last Updated: Aug. 9, 2019

Greg Pintilie

gregp@slac.stanford.edu

Overview

- MapQ is a UCSF Chimera plugin that allows calculation and visualization of Q-scores, which measure resolvability of atomic features in cryoEM maps using a fitted model.
- This tutorial covers the following topics:
 1. Get the software (UCSF Chimera and the MapQ plugin)
 2. Getting an example map and model
 - Feel free to use your own map and model here
 3. Use the MapQ interface to calculate and visualize Q-scores for amino acid residues and solvent atoms.
- The plugin has been tested on MacOS (10.13.16), Linux (Ubuntu 14.04), and Windows 10 with Chimera 1.12, 1.13 and 1.14.


1. Getting the Software

- Download Chimera for your system from
 - <https://www.cgl.ucsf.edu/chimera/download.html>
 - Either production or daily release should work.
- Download and install the MapQ plugin
 - <https://github.com/gregdp/mapq>
 - Follow the Install instructions at the above link.


2. Loading the map and model

- Download the map and model from the following link:
 - https://github.com/gregdp/mapq/blob/master/data/QScore_Apoferritin_Tutorial.zip
 - Once downloaded, unzip the file. It should contain two files, a map and a model:
 - Emd20026_prot.mrc – this is a portion of the entire Apoferritin map representing a single protein and solvent atoms around it.
 - 3ajo_fitted_refined.pdb – this is the X-ray model PDB:3ajo which was first rigidly fitted to the cryoEM map, and then refined using Phenix.real_space_refine
- Open Chimera
- Load the map:


- Open ModelZ dialog:
 - Tools menu -> Volume Data -> MapQ
- The ModelZ Dialog looks like this:



- A map, model, and chain are already selected.
 - If not, or to select a new map/model/chain, click on the respective fields just under the sequence.
- The sequence for the selected chain is shown
 - The sequence is made up of amino acid residues (e.g. Tyrosine, T, Serine, S, Glutamine, Q, etc.)
- Residues that are in 'loops' are shown with gray background.
- Residues that are in 'helices' are shown with reddish background.
- Select an entire loop/helix (Ctrl+Click+Drag on the sequence).
 - This selects just that part of the protein, and,
 - The map around the selection will be extracted and shown.
- With each part selected, click 'At.', 'Rib.', 'SCs', '~SCs' to show/switch between seeing Atoms, Ribbon, Side Chains, No-Chains. This will show the protein in different ways, as illustrated below.




- Make sure '**Ligands**' is checked in the MapQ dialog (to the right of 'Sequence Select:'). Then, select the entire second helix. The selection and extracted map will now include nearby 'ligand' atoms, which in this case are water and Mg atoms as below:




- Note how at lower thresholds you can see contours around some of the Mg and water (O) atoms. These atoms were placed in the X-ray model, and are also resolved in the cryoEM map.

- To calculate Q-scores, you can press either:
 - “**Calc**” button
 - This will calculate Q-scores for all atoms in the selected model and chain. Select ‘All’ in the Chain field to calculate for all atoms in the model.
 - For this protein, it will take about 4 min to calculate the Q-scores for each atom. You can see the progress and ETA at the bottom of the main Chimera window. To cancel, you can press the X button.
 - “**Load**” button
 - Once the Q-scores are calculated once, a new file will appear in the folder of the map, which is the same as the original model, but will have Q-scores stored in the B-factor column.
 - Pressing this button will look for this file and load it instead of calculating Q-scores again.
 - “**Sel**” button
 - This will calculate Q-scores only for atoms that are selected (with Ctrl+Left Click – press the up arrow to expand selection to residue, secondary structure, or entire chain).
 - It typically takes just a few seconds to calculate the Q-score for each atom, so this is faster than calculating Q-scores for an entire model.
- After calculating Q-scores using the **Calc** button (or loading them with the **Load** button), you will see the area above the sequence show a bar graph of Q-scores for each residue.
 - On the bottom are backbone Q-scores, on the top are side chain Q-scores:



- Because this map is very good, you will see most residues in this case have good Q-scores – the bars for these are green.
 - A couple have lower (red). These tend to be residues which are more dynamic (e.g. Asp, Glu), or portions of loops.
- To visualize Q-scores on the model, you can press a button to the right of the Color: label:
 - **Backbone** – colors the model ribbon by backbone score (green is high, red is lower Q-score).
 - **Side Chain** – colors ribbon by side chain score.
 - **Atoms** – shows a label with the Q-score next to each atom, as below
 - **Random** – colors each chain a random color



- If you select the '**Keep**' checkbox, when selecting a new part of the sequence, the previous selection is also included in what is shown.
 - So for example, you can extract/show two parts of the model at the same time.
 - Press the '<' button to undo the last selection (helpful if the new selection was not what you were looking for).

- Using this mode, we can show, for example, the GLU 27 residue which is close to the Mg on the right, while keeping the previous selection:

